

To study the effect of concentration on the rate of reaction between sodium thiosulphate and hydrochloric acid. (b) To study the effect of temperature on the rate of reaction between sodium thiosulphate and hydrochloric acid.

Video Reference :- <https://www.youtube.com/watch?v=dGfJSuJW8F8>

Materials Required :-

- 1 Molar Sodium Thiosulphate Solution
- 1 Molar Hydrochloric Acid
- 5 250 ml Conical Flask (A, B, C, D, E)
- White Tile With X Mark In Centre
- Beaker (Distilled Water)
- Measuring Cylinder

Procedure :-

- Conical Flash (A) Should Contain 10 ml of 1 Molar Sodium Thiosulphate Solution.
- Conical Flash (B) Should Contain 20 ml of 1 Molar Sodium Thiosulphate Solution.
- Conical Flash (C) Should Contain 30 ml of 1 Molar Sodium Thiosulphate Solution.
- Conical Flash (D) Should Contain 40 ml of 1 Molar Sodium Thiosulphate Solution.
- Conical Flash (E) Should Contain 50 ml of 1 Molar Sodium Thiosulphate Solution.