

To prepare crystals of

a) Ferrous ammonium sulphate (Mohr salt), $\text{FeSO}_4 \cdot (\text{NH}_4)_2\text{SO}_4 \cdot 6\text{H}_2\text{O}$.

Video Reference :- <https://www.youtube.com/watch?v=Kqk7kC03ofw>

Materials Required :-

- 7 Gram Ferrous Sulphate
- 3.5 Gram Ammonium Sulphate
- Dilute Sulphuric Acid (Glass Bottle)
- Alcohol (Plastic Can)
- Distilled Water
- Cold Water
- China Dish
- Funnel
- Spatula
- Glass Rod
- Beakers
- Measuring Cylinder (10 ml, 100 ml)
- Crystallising Dish
- Trough
- Watch Glass
- Filter Paper
- Stand with Clamp
- Hot Plate

b) Potash alum $\text{K}_2\text{SO}_4 \cdot \text{Al}_2(\text{SO}_4)_3 \cdot 24\text{H}_2\text{O}$

Video Reference :- <https://www.youtube.com/watch?v=bhoMvPJKc24>

Materials Required :-

- 2.5 Gram Potassium Sulphate
- 10 Gram Ammonium Sulphate
- Dilute Sulphuric Acid (Glass Bottle)
- Alcohol + Cold Water (Plastic Can)
- Distilled Water
- Cold Water

- China Dish
- Spatula
- Glass Rod
- Beakers
- Measuring Cylinder (10 ml, 100 ml)
- Crystallising Dish
- Watch Glass
- Hot Plate

c) **Potassium ferric oxalate $K_3[Fe(C_2O_4)_3] \cdot 3H_2O$.**

Video Reference :- <https://www.youtube.com/watch?v=SilzJBCEins>

Materials Required :-

- 3.5 Gram Anhydrous Ferric Chloride
- 4 Gram Potassium Hydroxide
- 4 Gram Hydrated Oxalic Acid
- 5.5 Gram Hydrated Potassium Oxalate
- Ethyl Alcohol
- Distilled Water
- Cold Water
- China Dish
- Buchner Funnel
- Spatula
- Glass Rod
- Beakers
- Measuring Cylinder (10 ml, 100 ml)
- Trough
- Watch Glass
- Filter Papers
- Black Paper
- Sand Bath